

Quick Start Guide

EM Evidence Sim · for instructors running scenarios · updated 1 October 2026

What you need

You can run a simulation with one laptop and one screen the team can see. Everything runs in a web browser at wmebemsim.netlify.app; there is nothing to install.

Device	Who uses it	Needed?
Laptop or tablet: the instructor screen	You, out of the team's sight	Always
A second screen, TV, projector or device: the room monitor	The team	Always
A tablet: the defib	The team	Only for Defib Sim, or if you want the team to use a defibrillator

- Every device needs an internet connection. The screens link to each other online using a session code.
- Sound (alarms and beeps) plays on the room monitor by default, so turn its volume up.
- You do not need to sign in to run the built-in scenarios.

Step 1: Open the app

1. On your laptop, open wmebemsim.netlify.app.
2. Note the six-character **Session Code** (for example K7PQ3M). It is shown on the setup screen and in the top right corner as "Session ID".
3. The code is what links the room monitor and the defib to your screen. Each group can keep the same code all day.

Starting with a new group? Press **New code**. Any monitor or defib still on the old code will need to join again.

Step 2: Choose a scenario

Pick a tab on the setup screen. For most teaching, use **Premade**.

Tab	Use it for	How to start
Premade	A ready-written scenario (anaphylaxis, sepsis, trauma, paediatrics and more)	Pick a category, press Load , read the brief, press Start Scenario

Tab	Use it for	How to start
Quick Sim	A blank patient you control by hand, for ad-hoc bedside teaching	Set age, weight and starting rhythm, press Start Quick Sim
Defib Sim	Defibrillation, cardioversion and pacing skills on the defib tablet	Pick Education or Assessment, pick a scenario, press Start Defib Sim
Random	A generated patient from a category and acuity you choose	Press Generate Scenario , then start it
Custom / Builder	Scenarios you have written or imported	Load or build one, then start it

The brief shown before **Start Scenario** has the patient story, equipment, expected progression, key interventions and learning objectives. Read it before the team comes in.

Step 3: Set up the room monitor (and defib)

The setup screen starts by asking **How are you running the sim?** This computer is your controller. Pick how many patient screens (1, 2 or 3) and where they are; this device remembers it.

1 screen shows the monitor and switches to the defib or ventilator when you bring them in. **2 screens:** the monitor, and the defib and ventilator. **3 screens:** monitor, defib and ventilator each on their own.

Plugged into this computer: press **Open screen 1** (and 2, 3). Each opens in its own window; drag it onto its screen and make it full screen.

Separate devices: scan each screen's QR code with that device. No camera? On the device, open the site, choose **This device is a patient screen**, type the code, pick what it shows and press **Connect**.

Then:

1. On the room monitor, tap **Tap to Enable Sound**. Browsers block sound until someone taps the screen.
2. Check the badge at the top of your instructor screen. It should read **Live** and name the linked screen (for example "Live · patient monitor").

Step 4: Run the scenario

You play the patient. When the team does something, you press the matching button, and the monitor and patient respond.

When the team...	You press...
Puts on ECG, sats, cuff and temperature probe	Attach on the Monitoring & access bar. Until then the room monitor shows no numbers
Asks for a blood pressure	Cycle NIBP Now . The room monitor only shows a new BP when the cuff cycles. When you change the BP, choose Change, don't send (they see it at the next cycle) or Change and cycle cuff now
Gives a drug or does a procedure	Find it in the search box or the tabs (Common, Drugs, Airway...) and tap it
Does something that should change a vital sign	Tap that vital's tile, type the new value, choose Now , 30 s, 2 min or 5 min, then Confirm (BP starts on Now and has its own two buttons)
Lets the patient arrest	ARREST ▼ and choose the rhythm
Gets the patient back	ROSC ▼

Also useful:

- **PAUSE / START** (or the space bar) stops and starts the clock.
- **Patient condition** shows whether the patient gets worse on their own (AUTO) or only when you change things (MANUAL). **Better** and **Worse** nudge the whole patient.
- **Event log** (bottom right) records everything with a time. Type a note and press **Add note**, or tap the flag beside any entry to discuss it later.
- **Simple / Full** (top bar): Simple hides the advanced controls. Use Full when you need them.

Step 5: Finish and debrief

1. Press **Finish** and confirm. The debrief screen opens on your laptop.
2. Use it to lead the debrief:
 - **Performance Summary**: which learning objectives were met, partly done or not done.
 - **Obs, interventions and changes**: a graph of the patient's obs over the whole session, with the rhythm, every treatment, every change you made and your flags marked and numbered. Hover over it to read the values at any moment.
 - **Simulation Log**: every action with its time.
3. Press **Download Report** to keep a copy, or **Print. Certificate** (Defib Sim only) makes a certificate of completion for the learner.
4. Press **Exit to Menu** to set up the next scenario. The room monitor and defib stay linked.

Quick fixes

Problem	Fix
Room monitor shows no numbers	The team has not "attached" monitoring yet. Press Attach on the Monitoring & access bar
No sound	Tap Tap to Enable Sound on the room monitor and check its volume. Tools → Sound plays on chooses which screen plays it
Badge says Sync error	Check the internet connection on each device
Room monitor or defib not following your changes	Check it has the same session code as your laptop; join again if not
Instructor tab closed by mistake	Reopen the site and press Resume under "Resume Previous?"
Defib Sim log times stuck at 00:00	The clock has not started. Press START , or it starts at the learner's first press on the defib

Keyboard shortcuts: **Space** = pause or start; **F** = finish; **D** = drug calculator; **T** = timer alerts; **?** = list of shortcuts. If a shortcut does nothing, click an empty part of the screen first.

For everything else, see the Full Instructor Guide.